## Delivering the Industry's Most Consistent Respirable Dose<sup>1</sup>

The AEROECLIPSE<sup>®</sup> II Breath Actuated Nebulizer (BAN) produces aerosol only when you breathe in. As a result, it gets small particles of medication deep in your lungs.


## Faster delivery of medications to the right place!

The AEROECLIPSE<sup>®</sup> II BAN gets medication quickly into your lungs. Now you will get all the medication you need to help you feel better - only you'll feel better faster!


## How can the AEROECLIPSE<sup>®</sup> II BAN change your treatments?

The AEROECLIPSE<sup>®</sup> II BAN allows your health care provider many treatment options that can include shorter treatment times, less medication put in your nebulizer or running the AEROECLIPSE<sup>®</sup> II BAN continuously. The result is better care that allows you to take your treatment in less time.

# AEROECLIPSE<sup>®</sup> II BREATH ACTUATED NEBULIZER

AE<sup>®</sup> II BAN use shows medication gets deep into the lungs.<sup>2</sup>

Conventional nebulizer use shows large particles ending up in stomach, not lungs.

Gamma Camera Deposition Images. Sangwan S, et al, Journal of Aerosol Medicine 2001.

## Instructions for use

Close your lips completely over the mouthpiece and do not cover the valve on the bottom of the mouthpiece. If using a mask, fit the mask tightly against your face. Take slow, deep breaths; try to hold each breath for 3 to 5 seconds. The green button on the top of the AEROECLIPSE<sup>®</sup> II BAN will go down when you breathe in and will go up when you breathe out. Be sure to also follow any additional instructions your health care provider prescribes. Please refer to the AEROECLIPSE<sup>®</sup> II BAN package insert for complete instructions and care.

Words or phrases accompanied by <sup>™</sup> and <sup>®</sup> are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2006-2016 Monaghan Medical Corporation. <sup>1</sup> trademarks and registered trademarks of the respective companies.

P/N 76198-07, 09/2017

<sup>1</sup> Investigation of Medication Delivery from Small Volume Nebulizers (SVN) and a Breath-Actuated Nebulizer (BAN) Using *In-Vivo* Generated Breathing Profiles Judy Schloss<sup>1</sup>, Dominic P. Coppola<sup>2</sup>, Jason Suggett<sup>3</sup>, Vivian T. Wang<sup>3</sup>, Cathy Doyle<sup>3</sup> and Mark W. Nagel<sup>3</sup>. (1 University of MN Medical Center, Fairview, MN; 2 Monaghan Medical Corp., Syracuse, NY; 3 Trudell Medical International, London, Canada.)

<sup>2</sup> This study used the AeroEclipse<sup>®</sup> Breath Actuated Nebulizer not the AeroEclipse<sup>®</sup> II Breath Actuated Nebulizer which is currently on the market. The AeroEclipse<sup>®</sup> BAN and AeroEclipse<sup>®</sup> II BAN are equivalent in regards to *in vitro* performance. "Are first and second generation, mechanically-operated breath-actuated nebulizers (BAN) comparable based on *in vitro* performance?" Schmidt J, Pevler J, Doyle C, Wiersema K, Nagel M, Mitchell J. Respiratory Drug Delivery, 2006:817-819.

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