## STRIVE® PEAK FLOW METER (PFM)

**1** Reset the flow indicator by gently sliding it to bottom of the scale until it stops. Hold STRIVE® PFM with your right hand by the grooved handle, ensuring the flow indicator is not obstructed.



**2** Stand or sit up as straight as possible, while still being comfortable.

**3** Take as deep a breath as possible. Do not inhale through the device.

**4** Put the STRIVE® PFM mouthpiece into your mouth, make sure you keep your lips tightly sealed around the mouthpiece.

**5** Blow out as hard and as fast as possible.

**6** Repeat steps 2 through 5 as many times as instructed by your physician or health care professional.

### Using your daily record:

The daily record covers 35 days. Write your name on the line provided, start with today's date, enter dates of next 35 days across the top of the record in spaces provided.

Under each date is a graph for charting Peak Expiratory Flow (PEF) up to three times a day. Begin entering results in your Daily Record in the left hand column of the chart. Subsequent results are entered to right.

Once the red flow indicator comes to rest at your final PEF, make an "X" to record the measurement number under the appropriate Date and Time heading.

These are abbreviated instructions for use, please remember to read the complete instructions that are packaged with your device.

Words or phrases accompanied by ™ and ® are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2016 Monaghan Medical Corporation. † trademarks and registered trademarks of the respective companies.

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