

More Time to Breathe Easy

Increased aerosol suspension time gives everyone more time to breathe easy**



FLOWSIGnal® Whistle

to promote proper technique when transitioning from a mask to a mouthpiece device

Flow Dynamic Valve and Baffle System

prevents patients from exhaling into the chamber, ensuring the remaining medication is available for the next breath

Consumer Friendly

- Easier to read user instructions
- Visually appealing design



Anti-Static Chamber

provides consistent aerosol delivery right out of the package

Universal pMDI Adapter

fits commonly prescribed pMDIs for easy insertion and removal of pMDIs



More Time to Breathe Easy.




monaghan™

AeroChamber
Plus **Z STAT™**

Anti-Static Valved Holding Chamber

Combining over 25 years experience in aerosol delivery and chamber design



In addition to being in vitro tested with all commonly prescribed metered dose inhaler medication used for respiratory conditions, AeroChamber® VHC has also been extensively tested In Vivo with Flovent†, Alvesco†, QVAR†, Serevent†, Albuterol†**

Most recommended chamber by leading metered dose inhaler companies*

#1 Physician prescribed brand of chamber*

Most peer-reviewed studies worldwide*

Ordering Information	
Part #	Description
79710	AeroChamber Plus® Z STAT® aVHC, case 10
79750	AeroChamber Plus® Z STAT® aVHC, case 50
80710	AeroChamber Plus® Z STAT® aVHC w/LRG ComfortSeal® Mask, case 10

For more information or a demonstration, contact your Monaghan sales rep:



monaghan™

Monaghan Medical Corporation, 5 Latour Ave., Suite 1600,
Plattsburgh, NY 12901, Customer Service 800-833-9653
www.monaghanmed.com

™ and ® are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation © 2010 Monaghan Medical Corporation
† trademarks and registered trademarks of the respective companies
*studies on file; ** AeroChamber® In Vivo Summary, PN76374

AeroChamber
Plus **Z STAT™**

Anti-Static Valved Holding Chamber