



Establishes target zones

The Monaghan Advantage

- ❖ Convenient, user-friendly device for reliable airway clearance
- ❖ Suitable for thoracic post-operative recovery; prevention or treatment of atelectasis; and lung volume expansion therapy
- ❖ Develops proper technique and actively engages patient in therapy
- ❖ Patient receives immediate feedback with a visual goal
- ❖ Achieve therapeutic pressure every time to maximize treatment effectiveness

MANOMETER ADAPTER FOR THE AEROBIKA® (OPEP)

The AEROBIKA® Oscillating Positive Expiratory Pressure (OPEP) Therapy System with the manometer adapter provides visual feedback during a patient's exhalation. The AEROBIKA® OPEP is convenient, user-friendly and provides clinicians and their patients with a means for reliable, drug-free, natural airway clearance.

The AEROBIKA® OPEP system opens weak or collapsed airways by creating positive pressure and oscillations simultaneously, which helps move mucus to the upper airways where it can be coughed out.

PROVIDING VISUAL FEEDBACK

THE FACTS

Proprietary pressure-oscillation dynamic delivers intermittent resistance and creates positive pressure and oscillations simultaneously. This opens weak or collapsed airways to mobilize and assist mucociliary clearance to the upper airways, where it can be coughed out.

Resistance is not position dependent, allowing the user to hold the device in the most natural position.

ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
62510	AEROBIKA® Oscillating Positive Expiratory Pressure (OPEP) Therapy System	cs/10
62610	AEROBIKA® Oscillating Positive Expiratory Pressure (OPEP) Therapy System w/Manometer	cs/10

Contact your Monaghan Sales Territory Manager for more information or a demonstration.



REGISTERED
ISO 13485
QUALITY SYSTEM

