

**What makes the
AEROECLIPSE® XL
Breath Actuated Nebulizer
(BAN) different?**

With breath actuated delivery, aerosol is only produced when you breathe in. This means medication is not wasted between breaths or if you need to take a break (for example if you need to cough or have a conversation). This puts you in control of your aerosol treatment and creates a safer home environment. Other nebulizers continuously produce aerosol regardless of whether you are inhaling, exhaling or resting, resulting in medication being lost to the air instead of delivered to your lungs.



**AEROECLIPSE® XL BAN
AND OMBRA® COMPRESSOR**



What is a Nebulizer?

A nebulizer is a medical device designed to work with a compressed air source to convert liquid medication into a fine mist so it can be inhaled directly into the lungs. You may need nebulizer treatments if you have asthma or lung diseases, such as cystic fibrosis or COPD.

If you have questions about the prescribed medication inside the nebulizer, please talk to your health care provider.

