



## VersaPAP® DEVICE

### Positive Airway Pressure Therapy (PAP)

PAP devices use pressurized air to help maintain positive airway pressure as the patient inhales and exhales through the device.

#### PAP Therapy And Outcomes<sup>1</sup>



Increased amounts of sputum produced



Improved vital signs (blood gas values, oxygen saturation, respiratory rate or pulse)



Improvement in breath sounds



X-ray indication improvement in atelectasis or gas trapping



Patients reporting ease in bringing up mucus and greater volumes of mucus

# PAP Therapy Is Helpful Because<sup>1</sup>...

Reduces Air Trapping

Prevents Or Reverses Atelectasis

Helps With Mucociliary Clearance

Improves Delivery Of Inhaled Medications



## ON INHALATION


- Increases the amount of air being inhaled
- Decreases the amount of work it takes to inhale
- Positive pressure helps stent airways open



## ON EXHALATION

- Takes more effort against the flow of air
- Keeps the smaller airways from collapsing before the next breath in

## What Is The Pressure Generated At Different Flows?<sup>2</sup>



## ORDERING INFORMATION

PART #	NAME/DESCRIPTION	UNITS PER CASE
68510	VersaPAP® Positive Airway Pressure Device with EZ TWIST Tubing	cs/10
68500	VersaPAP® Positive Airway Pressure Device with EZ TWIST Tubing	cs/100
68310	VersaPAP® Positive Airway Pressure Device with Universal Tubing	cs/10
68300	VersaPAP® Positive Airway Pressure Device with Universal Tubing	cs/100

Contact your Monaghan Sales Territory Manager for more information or a demonstration.

518.561.7330 phone Monaghan Medical Corporation  
 800.833.9653 toll free Plattsburgh, New York 12901  
 518.561.5088 fax [monaghanmed.com](http://monaghanmed.com)

monaghan means it matters®

References: 1. AARC Clinical Practice Guideline: Use of Positive Airway Pressure Adjuncts to Bronchial Hygiene Therapy. Respir Care 1993;38(5):516-521. 2. Coppolo D, et al. Am J Respir Crit Care Med 2019;199:A4696.

Words or phrases accompanied by™ and \* are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2019 Monaghan Medical Corporation. 68503-01, 12/2019