



**AEROCHAMBER2GO® VHC
RESCUE/RELIVER INHALER
COMPATIBILITY LIST**

Explore the current list of Rescue/Reliver inhalers that can be used and stored within the AEROCHAMBER2GO® VHC

**Aero2go®
Chamber 2go**

Anti-Static Valved Holding Chamber
For Use With Metered Dose Inhalers

RESCUE INHALERS IN THE US (INCLUDING SMART/MART)

<div>INHALERS STORED IN UPRIGHT POSITION</div> <div></div>	<div>SABA</div> <div></div> <div>VENTOLIN[†]</div> <div>GSK</div>	<div>SABA</div> <div></div> <div>ALBUTEROL[†]</div> <div>PRASCO</div>	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>LUPIN</div>	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>ARMSTRONG/AMPHASTAR</div>
	<div>SABA</div> <div></div> <div>ATROVENT[†]</div> <div>BOEHRINGER INGELHEIM</div>	<div>SMART/MART</div> <div></div> <div>DULERA[†]</div> <div>MERCK</div>	<div>SMART/MART</div> <div></div> <div>BREYNAT[†]</div> <div>VIATRIS</div>	<div>OTC</div> <div></div> <div>PRIMATENE MIST[†]</div> <div>ARMSTRONG</div>
<div>INHALERS STORED IN INVERTED POSITION</div> <div></div>	<div>SABA</div> <div></div> <div>PROVENTIL[†]</div> <div>MERCK</div>	<div>SABA</div> <div></div> <div>PROAIR[†]</div> <div>TEVA</div>		
	<div>SMART/MART</div> <div></div> <div>AIRSUPRA[†]</div> <div>ASTRAZENECA</div>	<div>SMART/MART</div> <div></div> <div>SYMBICORT[†]</div> <div>ASTRAZENECA</div>	<div>SMART/MART</div> <div></div> <div>GENERIC SYMBICORT[†]</div> <div>ASTRAZENECA</div>	
<div>INHALERS CAN BE USED WITH CHAMBER</div> <div></div> <div>BUT NOT STORED INSIDE</div>	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>CIPLA</div>	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>SANDOZ</div>	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>HIKMA</div>	<div>SABA</div> <div></div> <div>XOPENEX[†]</div> <div>LUPIN / TEVA</div>
	<div>SABA</div> <div></div> <div>ALBUTEROL</div> <div>PADAGIS</div>			

The "Rescue/Relievers" shown above are representative of the MDIs known to be available in the USA as of December 2025, and are not inclusive of all inhalers compatible with the AEROCHAMBER2GO® VHC. Words or phrases accompanied by "™" and "®" are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2025 Monaghan Medical Corporation. [†] trademarks and registered trademarks of the respective companies.