
















AEROCHAMBER2GO® VHC RESCUE/RELIEVER INHALER COMPATIBILITY LIST

Explore the current list of rescue/reliever inhalers that can be used and stored within the AEROCHAMBER2GO® VHC

Aero2go™ Chamber

Anti-Static Valved Holding Chamber
For Use With Metered Dose Inhalers

RESCUE INHALERS IN THE US (INCLUDING SMART/MART)

INHALERS STORED IN UPRIGHT POSITION 				
	VENTOLIN[†] GSK	ALBUTEROL[†] PRASCO	ALBUTEROL LUPIN	ALBUTEROL ARMSTRONG/AMPHASTAR
				
	ATROVENT[†] BOEHRINGER INGELHEIM	DULERA[†] MERCK	BREYNAT[†] VIATRIS	PRIMATENE MIST[†] ARMSTRONG
INHALERS STORED IN INVERTED POSITION 				
	PROVENTIL[†] MERCK	PROAIR[†] TEVA		
				
	AIRSUPRA[†] ASTRAZENECA	SYMBICORT[†] ASTRAZENECA	GENERIC SYMBICORT[†] ASTRAZENECA	
INHALERS CAN BE USED WITH CHAMBER 				
	ALBUTEROL CIPLA	ALBUTEROL SANDOZ	ALBUTEROL HIKMA	XOPENEX[†] LUPIN / TEVA
				
	ALBUTEROL PADAGIS	ALBUTEROL AMNEAL		

The "Rescue/Relievers" shown above are representative of the MDIs known to be available in the USA as of July 2026, and are not inclusive of all inhalers compatible with the AEROCHAMBER2GO® VHC. Words or phrases accompanied by "™" and "®" are trademarks and registered trademarks of Monaghan Medical Corporation or an affiliate of Monaghan Medical Corporation. © 2025-2026 Monaghan Medical Corporation. "†" trademarks and registered trademarks of the respective companies. "™" trademarks and registered trademarks of the respective companies.